

Pedigo Hydraulic Surgeon's Stool | P-7000



Stock Photo - Options May Vary

FACTORY NEW

The Pedigo P-7000 Hydraulic Surgeon's Stool is engineered for use in surgery, anesthesiology, ENT, and reconstructive procedures. Constructed using Pedigo's STERIL-GARD philosophy, the stool minimizes areas where dirt, fluid, and bacterial debris can become entrapped, reducing bacterial transmission.

The P-7000 features foot-operated hydraulic height adjustment with a single pedal that ranges from 20 3/4" to 28 1/4". The specially contoured Sensafoam seat and back provide ergonomic support during extended procedures. A chrome-plated 6-caster base with three locking casters ensures stability and mobility.

All components are designed for easy cleaning and long-term durability. The stool includes height-adjustable and articulating arm/wrist-rests, and the seat rotates 360 degrees independently of height adjustment.

FEATURES

- Foot-Controlled Hydraulic Operation
- Single Pedal Height Adjustment
- Steril-Gard Design
- Ideal for Surgery, Ophthalmology, Anesthesiology, ENT, and Reconstructive or Plastic Surgery
- Specially Contoured Seat for Added Comfort
- Height Adjustable and Articulating Arm/Wrist-Rests
- Chrome-Plated 6-Caster Base
- 3 Locking Casters

SPECIFICATIONS

- Foot-Controlled Hydraulic Operation
- Base Diameter: 21"
- Seat Height Adjustment Range: 20 3/4" - 28 1/4"
- Seat Rotates 360° Independently of Height
- Casters: 3"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)