

# Pedigo Infusion Pump Stand for ivNow Thermal Fluid Warmer | P-1080-6



Stock Photo - Options May Vary

## FACTORY NEW

The Pedigo IP-1080-6 ivNow Thermal Fluid Warmer Infusion Pump Stand is designed specifically for ambulatory infusion therapy applications. The P-1080-6 features accommodate up to six infusion pumps, allowing for efficient use in high-demand environments.

Additionally, a Bag Drop Prevention Device ensures the controlled descent of the pole, preventing accidental drops. The stand includes a durable stainless steel pole with hooks for added convenience and durability.

The IP-1080-6 also features heavy-duty 3" sealed casters, ensuring easy mobility and stability during use in any clinical environment. These thoughtful design features, combined with its sturdy construction, contribute to the stand's overall functionality.

## FEATURES

- Designed specifically for ambulatory infusion therapy applications
- Accommodates up to six infusion pumps
- Bag Drop Prevention Device allows for controlled descent of pole
- Stainless steel pole and hooks
- Adjustable height is the highest in the industry!

## SPECIFICATIONS

- Accommodates up to six infusion pumps
- Construction: Stainless steel pole, base, and hooks
- Adjustable height: from 69" to 124"
- Heavy-duty 3" sealed casters
- Finish: Extremely durable powder coat finish

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)