

# Pedigo iStand Infusion Pump Stand w/ 6 Hook | P-1080



Stock Photo - Options May Vary

## FACTORY NEW

Pedigo P-1080 iStand Infusion Pump Stand w/ 6 Hook. Pedigo designed the iStand specifically for ambulatory infusion therapy applications and accommodates up to six infusion pumps. The pole uses stainless steel that resists staining, corrosion, and rusting. Features Brass Bushing.

It provides smooth operation while reducing galling and wear. Features: The base features a low center of gravity design, and the 27 lbs. base offers improved stability.

The iStand's small footprint of only 23" reduces stumbles and falls while the patient ambulates. The base has an Extremely Durable Powder Coat Finish. Pedigo's NEW Infusion Pump Hoop (see options below) measures 29"H x 19" W and allows multiple pumps to be installed.

## FEATURES

- The stainless steel pole and hooks resist staining, corrosion, and rusting
- Brass Bushing provides smooth operation while reducing galling and wear
- The NEW Infusion Pump Hoop measures 29"H x 19"W and allows for multiple pumps to be installed
- Hoop allows for rotational and vertical adjustment
- The base has a low center of gravity design and 27 lbs. base offers improved stability

## SPECIFICATIONS

- Adjustable Height: 69" to 124"
- Casters: Heavy-duty, 3" Sealed Casters
- Base Finish: Durable powder coat
- Construction: Stainless Steel
- Base Weight: 27 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)