

Pedigo Mayo Stand/Medical Tray Unit | P-65



Stock Photo - Options May Vary

FACTORY NEW

The Pedigo Mayo Stand/Medical Tray Unit Model P-65 features a durable all-welded stainless steel construction with a chrome-plated frame, ensuring long-lasting performance and a sleek, professional appearance. Adjust the height (30" to 40") using the hand-operated friction lock for precise positioning during procedures.

The unit features 2" double ball bearing casters for smooth mobility and stability. The removable stainless steel tray provides a convenient surface for organizing surgical tools and supplies.

This versatile stand ensures optimal functionality and ease of use in busy medical environments. Note: Click the following link to view the wide variety of Mayo Stand models available to best meet your needs from Venture Medical.

FEATURES

- Stainless Steel Construction
- Chrome Plated
- Adjustable Hand-Operated Friction Lock
- Removable Tray

SPECIFICATIONS

- Construction: All-Welded Stainless Steel
- Removable Stainless Steel Tray Dimensions: 19-1/8" x 12-5/8"
- Hand Friction Lock Height Range: 30" - 40"
- Weight Capacity: 40 lbs. (evenly distributed)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)