

# Pedigo Narrow Transport/Trauma Stretcher | 7500-N



Stock Photo - Options May Vary

## FACTORY NEW

The Pedigo 7500-N Narrow Transport/Trauma Stretcher is designed for facilities where space is a consideration, with overall dimensions of 29.5"W x 83"L. It features a 750-pound load capacity and superior stability with Instant Steer Six Wheel Steering for improved maneuverability and control. The stretcher accommodates multiple clinical settings including OR, ED, and PACU.

The 7500-N includes full-length tuck-away side rails that create a zero transfer gap for easier patient transfer. The rail system wipes clean and disinfects easily. An ergonomic Quick-Release O2 Tank Holder mounts on the side and sits under the patient surface, reducing caregiver bending.

A sturdy steel center column with a single hydraulic lift cylinder uses sixteen precision ball bearings for smooth patient lift and lower. The stretcher offers Trendelenburg positioning, pneumatic assist Fowler backrest, and four corner IV pole receptacles for comprehensive patient care support.

## FEATURES

- 750 lbs. weight capacity
- Instant Steer 6 Wheel Steering (IS6)
- Side rails - the full-length, dual release
- Center column design with precision ball bearings
- Quick-Release O2 Tank Holder (US Patent 7,865,984 B2)
- Positive/negative Trendelenburg dual pneumatic assist cylinder
- Pneumatic assist Fowler (600N 135# force)
- Four corner IV pole receptacles

## SPECIFICATIONS

- Weight Capacity: 750 lbs
- Overall length: 83"
- Overall Width: 29.5"
- Height Range - Low: 23"
- Height Range - High: 32.75"
- Fowler Backrest: 0°-90°
- Trendelenburg: +12°
- Reverse Trendelenburg: -7°

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)