

Hill-Rom P8040 Trauma Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Hill-Rom P8040 Trauma Stretcher is engineered for high-pressure clinical environments, offering durable construction with smooth mobility, precise positioning, and advanced safety features. It supports up to 700 lbs and includes a full-length X-ray platform, IV pole, and Steering Plus™ system for enhanced control and maneuverability.

The stretcher features dual-locking brakes, ergonomic handles, and a radiolucent deck to ensure safe handling and quality imaging. Its power-washable frame is easy to clean and requires minimal maintenance, making it ideal for demanding clinical settings.

Trusted in trauma and emergency departments, the P8040 improves workflow and supports high-quality patient care in fast-paced environments.

FEATURES

- 700 lbs weight capacity
- Full-length X-ray lift platform accepts all standard-size X-ray cassettes and grids
- Next generation Steering Plus™ system for crisp cornering and control
- Dual-locking four-wheel brakes for a strong, stable platform
- Power washable
- Integrated IV pole transport device
- Ergonomic push handles
- Radiolucent deck

SPECIFICATIONS

- siderails stored: 29.25"
- siderails up: 32"
- height: 11"
- length: 47"
- 26-inch width

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)