

# Pedigo Wide Transport/Trauma Stretcher | 7500-W



Stock Photo - Options May Vary

## FACTORY NEW

The Pedigo 7500-W Wide Transport/Trauma Stretcher is designed for patient transfers across multiple clinical departments including OR, ED, PACU, and Radiology. It accommodates patients up to 750 pounds and features a center column design with precision ball bearings for smooth hydraulic operation. The full-length, dual-release side rails create a zero transfer gap for easier patient movement.

The stretcher incorporates Instant Steer Six Wheel Steering (IS6) with two wheels positioned directly under the chassis for improved maneuverability and stability. An ergonomic Quick-Release O2 Tank Holder mounts on the side beneath the patient surface, reducing caregiver bending. The rail system is designed for easy cleaning and disinfection.

Key positioning features include pneumatic assist Fowler backrest (0-90 degrees), positive and negative Trendelenburg capability, and a mechanical level lock. The stretcher includes four corner IV pole receptacles, four restraint strap positions, and 8-inch easy mobility casters with foot pedal-activated braking and steering.

## FEATURES

- 750 lbs. weight capacity
- Instant Steer 6 Wheel Steering (IS6)
- Side rails - full-length, dual release
- Quick-Release O2 Tank Holder (US Patent 7,865,984 B2)
- Positive/negative Trendelenburg
- Pneumatic assist Fowler (600 N 135# force)
- Four corner IV pole receptacles
- 8" easy mobility casters

## SPECIFICATIONS

- Weight Capacity: 750 lbs
- Overall length: 83"
- Overall Width: 33.5"
- Height Range - Low: 23"
- Height Range - High: 32.75"
- Fowler Backrest: 0°-90°
- Trendelenburg: +12°
- Casters: 8"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)