

Stryker 1020 Trauma Stretcher - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 1020 Trauma Stretcher is designed to meet the urgent demands of trauma care environments, offering strength, stability, and easy access for checking and treating seriously injured or immobile patients. Its design allows full-length X-rays without moving the patient, helping save time and reduce risk.

The stretcher features a secure side C-spine cassette holder and a slide-out tray for easy imaging access and clear results. A 90-degree backrest allows better positioning and supports patient care. Its strong build and trauma-ready features make it ideal for emergency rooms, trauma areas, and critical care transport.

Each unit is refurbished and tested by technicians to meet high-performance standards, ensuring reliable service and safety in critical settings.

FEATURES

- 4" New Mattress
- 26" wide patient surface
- Dual cylinder pneumatic backrest
- Fold-down siderails
- Four IV receptacles
- The full-length X-ray cassette tray
- Integrated storage tray with an oxygen bottle holder
- Four-wheel, steel-ring brakes with dual end activators

SPECIFICATIONS

- Model Number: 1020
- Overall Length: 83"
- Overall Width: 31.5"
- Weight Capacity: 500 lbs
- Height Range: (to litter top): High 38.5"
- Height Range: Low 24.75"
- Patient Surface: 26"W x 76"L
- Caster Diameter: 8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)