

Pedigo Over-Operating Specialty Table | P-5182



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Over-Operating Specialty Table model number P-5182. Pedigo's Over-Operating Specialty Table is stainless steel and features easy one-hand height adjustment.

It also features Total lock casters, brakes, swivel locks, and adjustable height from 37" to 52" high from the floor with 27¼" clearance between legs. Note: Click the following link to view our Complete Catalog containing a wide variety of Instrument Assembly Tables.

Looking for other Stainless Steel Operating Room Equipment to meet your specific needs? Click the link to view our Complete Catalog of Stainless Steel Equipment.

FEATURES

- Stainless steel
- Easy one-hand height adjustment
- Total lock casters, with brakes and swivel locks
- Adjustable height range from 37" to 52" high from the floor
- Additional Sizes and Configurations are available upon request

SPECIFICATIONS

- Casters: 3" Total Lock w/brakes and swivel locks
- Adjustable Height: 37" to 52" from floor
- The clearance between legs: 27 ¼"
- Overall Size: 34" W x 42" D

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)