

Pedigo Pediatric Stretcher/Crib Special Package | 500-SPEC



Stock Photo - Options May Vary

FACTORY NEW

The Pedigo Pediatric Stretcher/Crib 500-SPEC is a full-function transport crib designed for pediatric patients with a unique 2-post support system that provides 97% access to the patient. The crib features counterbalanced side rails with 2-inch bar spacing and a one-handed dual-action release handle that is childproof and easy to operate. Casters with central locks and steering controls from two corners enable increased maneuverability.

Key operational features include Trendelenburg and Reverse Trendelenburg positioning, full-function Fowler capability, and two IV pole receptacles. The side rails are counterbalanced to prevent sudden drops when released, protecting both child and provider from injury. The epoxy finish on the rails is warmer to the touch than chrome, and rounded corners enhance safety.

This special package includes the stretcher crib, a standard 3.5-inch mattress, utility shelf, and IV pole, providing a complete transport solution for pediatric care environments.

FEATURES

- Casters with a central lock and steering with controls from two corners
- Trendelenburg/Reverse Trendelenburg
- Two IV pole receptacles
- One-handed rail operation
- Rails with 2" bar spacing
- Full function Fowler
- Unique rounded corners
- Security Handle: Squeeze the Handle Together, Lift the Rail Upward or Downward

SPECIFICATIONS

- Casters: 5" with a central lock and steering with controls from two corners
- Finish: White powder coat
- Trendelenburg/Reverse Trendelenburg
- Rail Spacing: 2" bar spacing

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)