

# Pedigo Revolution™ Case Cart | RCC-242



Stock Photo - Options May Vary

## FACTORY NEW

Pedigo Revolution™ Case Cart: RCC-242. The Revolution™ closed case cart features an ergonomic design that makes it lightweight and easy to handle. New re•volve™ casters improve maneuverability and reduce push/pull force by 32% compared to standard case carts.

The doors feature clear, impact- and scratch-resistant poly windows. The aluminum doors come powder-coated with a custom high-temp, chip-resistant resin. The Re•volve™ Caster is made from glass-reinforced engineered thermoplastic.

Rated at 280° Fahrenheit, this high-temp caster features four precision bearings in each caster – the best in the industry. Twin precision sealed ball bearings improve cart handling and maneuverability and reduce the overall push/pull force by 32%.

## FEATURES

- Maintenance Free
- Caster is made from a glass-reinforced engineered thermoplastic
- Caster rated to 280° Fahrenheit
- Caster has four precision bearings in each caster and direction lock
- Caster improves maneuverability and reduces the overall push/pull force by 32%

## SPECIFICATIONS

- Outside Dimensions: 27½" D x 41½" H x 45½" W (top) / 46½" W (bottom)
- Inside Dimensions: 237/8" D x 26¾" H x 37¼" W
- Construction: Combination of Aluminum and Stainless Steel
- Windows: Scratch Resistant Poly
- Rated to 280 degrees

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)