

Pedigo Sealed Surgical Case Cart, Double Door | SCC-245



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Sealed Surgical Case Cart with Double Door model number SCC-245. Pedigo is proud to introduce Sealed Surgical Case Carts. These carts provide added protection against dust and debris, featuring a complete seal around the door perimeter and center stop.

Withstands cart wash/dry temperatures and sanitizing detergents. All Carts be built to your specifications.

Additional Features: All Surgical Case carts are built to Steril-Gard specifications, Fully-welded all stainless steel construction, Sound-proofed tops and doors for quiet operation, Exterior push handles for better control and more interior space, and a one-piece non-marring continuous wrap-around bumper. Note: Sealed Surgical Case Cart Shown with optional (sold separately) Accessories.

FEATURES

- Provide an added level of protection against dust and debris
- Full seal around the door perimeter and center stop
- Withstands cart wash/dry temperatures and sanitizing detergents
- All Surgical Case carts are built to Steril-Gard specifications
- Fully welded all stainless steel construction

SPECIFICATIONS

- Construction: All Stainless Steel
- 6" Chrome Plated Swivel Casters
- Overall Dimensions: 28-5/8"W x 28-5/8"D x 56"H
- Cabinet Dimensions: 23"□W x 24¼"□D x 47Ē"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)