

# Pedigo Specialty Back Table, Space Station | CDS-3060



Stock Photo - Options May Vary

## FACTORY NEW

Pedigo Specialty Back Table: 60" Wide Space Station CDS-3060. A specialty back table provides additional room for large cases such as orthopedics, spinal fusions, neuro, and craniotomies.

The tuck-away cantilevered shelf holds multiple trays and is angled for clear instrument observation. Users can easily roll the CDS-3060 into position using heavy-duty ball-bearing brake/swivel casters.

Additional CDS-3060 Features: Easy to lock this Specialty Back Table in place - engage the single-pedal brake and swivel caster locks. Strong, heavy-duty, all stainless steel construction. Pneumatic assist enables one-person operation in raising or lowering the cantilevered shelf. The shelf raises easily when extra space is required and lowers into its storage position when not needed.

## FEATURES

- Specialty back table designed with additional room
- For large cases such as orthopedics, spinal fusions, neuro, and craniotomies
- Tuck-away cantilevered shelf holds multiple trays and is angled for clear observation of instruments
- Easily rolled into position
- Heavy-duty ball-bearing brake/swivel casters

## SPECIFICATIONS

- Construction: Stainless Steel
- Main Table Dimensions: 30"D x 60"W x 34"H
- Top Shelf Dimensions: 20" D x 60" W
- Top Cantilevered Shelf Weight Capacity: 160 lbs
- Table Top Weight Capacity: 230 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)