

Pedigo Specialty Back Table, Space Station | CDS-3072



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Specialty 72" Back Table: Space Station CDS-3072. A specialty back table provides additional room for large cases such as orthopedics, spinal fusions, neuro, and craniotomies.

The tucked-away cantilevered shelf holds multiple trays and angles to clearly display instruments. Furthermore, the Space Station rolls easily into position with heavy-duty ball-bearing brake/swivel casters.

Additional CDS-3072 Features: Easy to lock this Specialty Back Table in place - simply engage the single-pedal brake and swivel caster locks.

FEATURES

- Specialty back table designed with additional room
- For large cases such as orthopedics, spinal fusions, neuro, and craniotomies
- Tuck-away cantilevered shelf holds multiple trays and is angled for clear observation of instruments
- Easily rolled into position
- Heavy-duty ball-bearing brake/swivel casters

SPECIFICATIONS

- Construction: Stainless Steel
- Main Table Dimensions: 30" D x 72" W x 34" H
- Top Shelf Dimensions: 20" D x 72" W
- Top Cantilevered Shelf Weight Capacity: 160 lbs
- Table Top Weight Capacity: 230 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)