

Pedigo SS Anesthetist Cabinet w/4 Drawers | P-1105-SS



Stock Photo - Options May Vary

FACTORY NEW

Pedigo Stainless Steel Anesthetist Cabinet w/4 Drawers, model number P-1105-SS. Pedigo offers a variety of anesthesia carts, cabinets, and treatment cabinets to fill the needs of many departments. Every model has smooth-rolling casters for easy mobility.

All stainless steel fully welded construction. The P-1105-SS is a Versatile Anesthetist cabinet that provides extra storage space and fully extending drawer hangers.

Additional Features Include All welded & polished stainless steel construction, Removable 12-5/8" x 19-1/8" stainless steel tray, and 3" double ball bearing swivel casters.

FEATURES

- Versatile Anesthetist cabinet provides extra storage space
- Full extending drawer hangers
- Removable 12-5/8" x 19-1/8" stainless steel tray
- All welded & polished stainless steel construction
- 3" double ball bearing swivel casters

SPECIFICATIONS

- Overall Dimensions: 19-1/8"W x 15-5/8"D x 33-3/8"H
- Three Drawers Dimensions: 13-1/2"W x 12"D x 3-1/2"H
- Lower Drawer Dimension: 13-1/2"W x 12"D x 7"H
- Rail Dimensions: 7/16 diameter x 1-1/2"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)