

Pedigo SS Anesthetist Table w/1 Drawer | P-1100-SS



Stock Photo - Options May Vary

FACTORY NEW

Pedigo SS Anesthetist Table w/1 Drawer: P-1100-SS. Pedigo offers a variety of anesthesia carts, cabinets, and treatment cabinets to fill the needs of many departments. Every model has smooth-rolling casters for easy mobility.

All stainless steel and fully welded construction. Ideal cart for storage of gloves, instruments, sutures, etc. Additional Features Include All welded and polished stainless steel construction. 2" double ball bearing swivel casters and 2 Removable Stainless Steel Trays, Towel Bar, and 2 Mask Hooks.

FEATURES

- Ideal cart for storage of gloves, instruments, sutures, etc.
- Two removable stainless steel trays
- Towel bar and two mask hooks
- All welded and polished stainless steel construction
- 2" double ball bearing swivel casters

SPECIFICATIONS

- Overall Dimensions: 21-3/8"W x 14"D x 34-1/4"H
- Trays: 2
- Tray Dimensions: 19-1/8"W x 12-5/8"D x 5/8"H
- Drawer Dimensions: 12-1/2"W x 15-3/8"D x 3-1/4"H
- Construction: Stainless Steel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)