

Pedigo SS Basin Stand w/Single Stainless Steel Basin | P-1078-SS



Stock Photo - Options May Vary

FACTORY NEW

The Pedigo SS Basin Stand w/Single Basin: P-1078-SS features a durable, rust-resistant stainless steel frame and basin, making it ideal for high-demand environments. The open frame design allows for easy access from all sides, providing convenience and efficiency.

This stand is equipped with smooth-rolling casters for easy mobility, allowing you to position it wherever needed. The stainless steel basin is easy to clean, offering a hygienic solution for medical, dental, or laboratory settings.

The combination of mobility, durability, and ease of maintenance makes this basin stand a reliable choice for any professional workspace.

FEATURES

- Stainless Steel Stand
- Stainless Steel Basin
- Open frame design for easy access from any side
- 2" diameter double ball bearing swivel casters
- Mobile

SPECIFICATIONS

- Basin Construction: Stainless Steel
- Stand Construction: Stainless Steel
- Overall Dimensions: 14D x 14"W x 34"H
- Basin: One removable 7-quart capacity stainless steel basin
- Basin Diameter: 13¼"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)