

Pedigo SS Double Basin Stand with Shelf | P-1079-W/S-SS

FACTORY NEW

The Pedigo SS Double Basin Stand with Shelf (Model P-1079-W/S-SS) features a durable stainless steel frame and double basins, offering a hygienic, long-lasting work surface. The open frame design provides easy access from all sides, and smooth-rolling casters ensure mobility.

The built-in shelf offers extra storage, while the stand's easy-to-clean design makes it ideal for medical, dental, or laboratory settings. Please call for pricing and additional details.



Stock Photo - Options May Vary

FEATURES

- Open frame design allows access from all sides
- Swivel casters
- Elevated ring basin bracket
- Mobile and easy to clean

SPECIFICATIONS

- Construction: Solid Stainless Steel
- Casters: 2" dia, double ball bearing, swivel casters
- Overall Size: 14"D x 28"W x 33 ¾"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)