

Pedigo Stackable SS Footstool Transport Cart | P-1015-TC



Stock Photo - Options May Vary

FACTORY NEW

The Pedigo P-1015-TC is a stainless steel transport cart designed for stackable footstools in high-traffic medical and clinical environments. Its all-welded steel construction with 1" tubular steel legs provides exceptional strength and durability while remaining tip-resistant and easy to clean with minimal maintenance requirements.

The cart accommodates up to six stackable footstools (model P-1015), offering a space-saving storage and transport solution. Swivel casters enable smooth mobility and easy maneuvering through tight or busy spaces, while the compact design keeps footstools readily accessible and organized.

Built for efficiency and practicality, the P-1015-TC simplifies footstool transport and organization in hospitals, clinics, and offices, maintaining a neat and efficient workspace.

FEATURES

- Strong, all-welded steel
- 1" tubular steel legs
- Tubular steel Handrail
- Tip resistant, rugged, and easy to maintain
- Swivel casters

SPECIFICATIONS

- 28"D x 14¾"W x 40½"H
- Weight Capacity: 300 lbs.
- Casters: 3" Swivel Casters
- 1" O.D. 16-gauge tubular steel
- Handrail Constructed: 1" O.D. tubular steel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)