

# Pedigo Stackable Stainless Steel Footstools | P-1015



*Stock Photo - Options May Vary*

## FACTORY NEW

Pedigo Stackable Stainless Steel Footstool model number P-1015: features strong, all-welded steel construction for durability. The footstool is tip-resistant, rugged, and easy to maintain.

Its legs are made from 16-gauge tubular steel for stability and strength. The corrugated black rubber top surface ensures a non-slip grip, while the rubber floor tips protect the floor and provide additional stability.

The stackable design allows for efficient storage when not in use.

## FEATURES

- Strong, all-welded steel
- Legs constructed of 16-gauge tubular steel
- Handrail constructed of tubular steel
- Tip resistant, rugged, and easy to maintain
- Corrugated black rubber top surface

## SPECIFICATIONS

- Overall Dimensions: 14"D x 19"W x 6"H
- Construction: All-welded steel
- Weight Capacity: 650 lbs.
- Leg Construction: 1" O.D. 16-gauge tubular steel

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)