

# Pedigo Stainless Steel Instrument Assembly Station | CDS-1000



Stock Photo - Options May Vary

## FACTORY NEW

Pedigo Stainless Steel Instrument Assembly Station: CDS-1000. Pedigo's Instrument Assembly Station is Designed to improve the efficiency of assembling cases in Sterile Processing—a wrap-around work surface for easy access to all areas. Features include electrically operated hydraulic height adjustment that provides comfortable, ergonomic height for all staff members.

Additional CDS-1000 Features: Three large pegboards for organizing instruments for storage and easy access and eight slide-out lower shelves for bins. Heavy-duty locking casters ensure easy mobility. Six removable stainless steel trays are recessed into the top work surface. A center drawer comes standard, with an optional keyboard tray.

Recessed trays hold small items like tags, labels, and container locks.

## FEATURES

- Designed to improve the efficiency of assembling cases in Sterile Processing. Wrap-around work surface for easy access to all areas
- Electrically operated hydraulic height adjustment provides comfortable, ergonomic height for all staff members
- Three large pegboards for organizing instruments for storage and easy access
- Eight slide-out lower shelves for bins
- 3" diameter, heavy-duty locking casters for easy mobility

## SPECIFICATIONS

- Casters: 3" diameter, heavy-duty locking
- Stainless steel construction
- Main Table Top Surface: 66" W x 30" D
- Height Adjustable: 36" to 47½" H

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)