

# Pedigo Stainless Steel 8.5 Quart Sponge Receptacle | P-1025-SS



Stock Photo - Options May Vary

## FACTORY NEW

The Pedigo Stainless Steel 8.5 Quart Sponge Receptacle: P-1025-SS offers a practical and durable solution for waste disposal in medical, dental, or laboratory settings. Designed with a balanced, tip-resistant structure, it ensures stability even during frequent use. The swivel casters provide easy mobility, allowing you to position the receptacle wherever needed.

The electro-polished stainless steel frame offers enhanced durability, corrosion resistance, and a sleek, professional appearance. With an 8 1/2 quart capacity, this receptacle is perfect for efficient storage and disposal of sponges and other materials. For additional details, refer to the tabs provided.

Note: Click the following link to view additional models, sizes, and features available, Kick Buckets.

## FEATURES

- Balanced, tip-resistant
- 2" swivel casters for easy mobility
- Electro-polished, stainless steel frame
- Perimeter bumper to protect walls
- Removable 8 1/2 Quart Basin

## SPECIFICATIONS

- Construction: Stainless Steel
- Basin Capacity: 8 1/2 Quarts
- Diameter: 16 1/4"
- Height: 8"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)