

Pedigo Surgeon's Stool with Procedure Rest | P-526

FACTORY NEW

The Pedigo Surgeon's Stool with Procedure Rest: P-526 is designed to provide comfort and functionality for medical professionals. The hydraulic system allows for smooth height adjustments, ensuring optimal positioning during procedures.

Full upholstery with special density foam offers lasting comfort and support, reducing fatigue during long surgeries or exams, while the durable construction ensures reliability in demanding environments. The stool rotates 360 degrees, offering enhanced mobility and flexibility for quick, efficient movement.

The procedure rest provides additional support, promoting a stable and ergonomic posture, which helps maintain comfort and efficiency throughout surgical and examination tasks. Ideal for clinical settings, this stool combines practicality with long-term durability.



Stock Photo - Options May Vary

FEATURES

- Weight capacity rating of 300 lbs
- Upholstered for comfort
- Rotates a full 360 degrees
- Special density foam for comfort and support
- Round, fully padded seat

SPECIFICATIONS

- Weight Capacity Rating: 300 lbs
- Rotation: Rotates 360°
- Seat Dimensions: 14½"
- Height Adjustment: 16¾" to 23¾"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)