

Pedigo Surgeon's Stool with Procedure Rest | P-6001

FACTORY NEW

The Pedigo Surgeon's Stool with Procedure Rest: P-6001 is designed to deliver comfort and mobility, making it an essential tool for medical professionals. The hydraulic lift allows for effortless height adjustments, ensuring the right positioning during various procedures.

Full upholstery and special density foam offer exceptional comfort and support, reducing fatigue even during long surgeries or exams. The stool rotates 360 degrees, allowing for flexible movement and quick adjustments without needing to reposition.

The procedure rest adds extra support, promoting a stable and ergonomic posture, which helps reduce strain and improve focus.



Stock Photo - Options May Vary

FEATURES

- Weight capacity rating of 300 lbs
- Rotates a full 360 degrees
- Upholstered for comfort
- Round, fully padded seat
- Height adjustment of 24-3/8" to 29½"

SPECIFICATIONS

- Seat Dimensions: 16"
- Height adjustment: 24-3/8" to 29½"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)