

# Pedigo Triple Pod ivNow® Modular Thermal Fluid Warmer | P-2303-A



Stock Photo - Options May Vary

## FACTORY NEW

The Pedigo P-2303-A Triple Pod ivNow® Thermal Fluid Warmer provides warmed fluid at your fingertips, exactly when and where you need it. These small, convenient units offer flexibility and you can place them in various combinations and environments.. Choose from countertop, wall-mounted, or mobile units based on your needs.

Features Include: ivNow quickly warms and maintains safe temperatures for intravenous and irrigation fluids. Additionally, the specially contoured warming surface cradles one (1) liter solution bag. Each cavity features individual control and an LED display for easy management.

Moreover, a sensor in the heating plate detects a bag's presence and engages the heating mechanism to warm the fluid.

## FEATURES

- Pedigo's ivNow fluid warmers are small, convenient units
- Can be placed in a number of different combinations and environments
- Countertop, wall mounted and mobile units are all available depending on your particular requirements
- Quickly warms and maintains safe temperatures of intravenous and irrigation fluids
- Each specially contoured warming surface cradles one (1) liter solution bag

## SPECIFICATIONS

- Dimensions: 18.1"H x 14.5"W
- Capacity: 3 bags (liter)
- Display: LED
- Temp Sensors: 2
- Accuracy: Within  $\pm 0.2^{\circ}\text{C}$  ( $\pm 0.3^{\circ}\text{F}$ ) of the setpoint temperature

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)