

Pedigo Wall Mounted ivNow®-4 Modular Thermal Fluid Warmer | P-2304-A



Stock Photo - Options May Vary

FACTORY NEW

Pedigo P-2304-A Wall Mounted ivNow®-4 Thermal Fluid Warmer. Pedigo's ivNow modular fluid warmers deliver warmed fluid at your fingertips, exactly when and where you need it. These small, convenient units mount on the wall for easy access.

Each unit holds four 1-liter bags and requires wall mounting using the provided brackets. Features Include: ivNow quickly warms and maintains safe intravenous and irrigation fluid temperatures.

Each specially contoured warming surface cradles one (1) liter solution bag. Additionally, each cavity has individual control and an LED display for precise temperature management.

FEATURES

- Pedigo's ivNow fluid warmers are small, convenient units
- Wall-mounted using Wall Mounting Bracket Kits
- Quickly warms and maintains safe temperatures of intravenous and irrigation fluids
- Each specially contoured warming surface cradles one (1) liter solution bag
- An individual control controls each cavity, and L.E.D. display

SPECIFICATIONS

- Dimensions: 24.2"H x 14.5"W Capacity: 3 bag (liter)
- Display: LED
- Temp Sensors: 4
- Accuracy: Within $+0/-2^{\circ}\text{C}$ ($+0/-3^{\circ}\text{F}$) of the setpoint temperature
- Note: Not approved for blood or blood products

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)