


Pedigo Wall Mounted ivNow®-5 Modular Thermal Fluid Warmer | P-2305-A



Stock Photo - Options May Vary

FACTORY NEW

Pedigo P-2305-A Wall Mounted ivNow®-5 Thermal Fluid Warmer offers reliable performance for medical environments. Pedigo's ivNow modular fluid warmers provide warmed fluid at your fingertips, exactly where and when you need it. Additionally, these small, convenient units mount on the wall.

To install, simply mount the five 1-liter bag capacity models using mounting brackets. Furthermore, they are designed for easy installation and quick access, ensuring your fluids stay at the right temperature for immediate use. As a result, you can count on fast, reliable fluid warming whenever necessary.

Immediate Availability at the Proper Temperature: ivNow quickly warms and maintains safe temperatures of intravenous and irrigation fluids.

FEATURES

- Pedigo's ivNow fluid warmers are small, convenient units
- Wall-mounted using Wall Mounting Bracket Kits
- Quickly warms and maintains safe temperatures of intravenous and irrigation fluids
- Each specially contoured warming surface cradles one (1) liter solution bag
- Each cavity is controlled by an individual control, and L.E.D. display

SPECIFICATIONS

- Dimensions: 29.6"H x 14.5"W Capacity: 5 bags (liter)
- Display: LED
- Temp Sensors: 2
- Accuracy: Within +0/-2°C (+0/-3°F) of the setpoint temperature
- Note: Not approved for blood or blood products

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)