

Replacement End Bumper Pads for Hausted 4160 Pediatric Crib



Stock Photo - Options May Vary

FACTORY NEW

Replacement end bumper pad set designed for the Hausted 4160 Pediatric Crib. This set includes two pads measuring 1" x 15.5" x 23" each, constructed with high-density foam for cushioning and support.

The pads are covered in durable Herculite Sure Check Comfort Fabric with welded corners to prevent fluid seepage. They feature flame and fluid resistance for enhanced safety in clinical environments.

The pads offer exceptional durability with resistance to bacteria, stains, and tears, making them suitable for frequent cleaning and long-term use in pediatric care settings.

FEATURES

- Fits End Bumper Pads for Hausted Pediatric Crib model number: 4160
- Pad thickness: 3 inches
- Set of two End Bumper Pads
- Upholstery Color: Black

SPECIFICATIONS

- (2) 1" x 15.5" x 23" Pads contain High Density foam
- Covered in durable Herculite Sure Check Comfort Fabric
- Welded corners to prevent fluid seepage
- Flame and Fluid resistance for safety

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)