

# Hausted SurgiChair ESC-EYE - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

The Hausted SurgiChair ESC-EYE is a multi-position chair and stretcher designed for patient treatment, transport, and recovery. The adjustable head section enables precise positioning, making it ideal for procedures requiring accuracy, particularly in ophthalmology, plastic surgery, and ENT applications.

This electric/battery-operated unit features an articulating headrest, padded armrests, and a fold-down adjustable footrest. The chair includes 2.5" cushions with a head pad, four-wheel brake and steer caster system, and a 24 V backup battery for continuous operation.

The package includes a hand-held remote control, new batteries, new pad set, and user manual. The thermoformed base cover and removable I.V. pole (optional) provide additional functionality for clinical environments.

## FEATURES

- New Pad Set
- 5" (127mm) casters
- Four-wheel brake and steer caster system
- 2.5" (64mm) cushions, including head pad
- Fold-down footrest (height adjustable)
- 24 VDC backup battery operation
- Electric/battery-powered handheld control
- Low profile articulating headpiece

## SPECIFICATIONS

- Height Range: High - 30.5" (768mm) to Low - 22.25" (569mm)
- Overall Width: 29.4" (747mm)
- Patient Surface: 22" (559mm)
- Overall Length: 75" (1905mm)
- ESC-EYE Model Input Voltage: 120 VAC, +5%, 60 Hz
- Amperage: Max. 3.5 A Duty Cycle: 10% six min. in 1 hr
- Padded Armrest Height: 12" (305mm)
- Padded Armrest Length: 28.375" (721mm)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)