

# Stryker 1089 Eye Surgery Stretcher – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker 1089 Eye Surgery Stretcher is specially designed for ophthalmic and head procedures, offering maximum head-end clearance to give surgeons the space needed for precise, delicate work. The stretcher features a locking work surface and mobile base to improve patient safety and streamline workflow during procedures.

The 1089 supports both superior and temporal surgical approaches with a dual-moving headpiece that provides better access during surgery. Built-in head controls tuck into the headpiece for easy movement and adjustment. A flexible tube system functions as both an air line and drape support, while dual-cylinder air backrest provides extra support during procedures.

The stretcher uses dual pedestal hydraulics for steady and reliable height changes, with controls positioned on both sides for smooth adjustment from either side. Additional features include fold-down siderails, four IV receptacles, four-wheel steel-ring brakes with foot-end brake and steer pedal, and a new 2-piece pad set.

## FEATURES

- New Patient Pad and Concave Headrest Pad
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual-side hydraulic controls
- Fold-down side rails
- Four IV receptacles
- Central-locking brakes with foot-end brake and steer pedal
- Dual-articulation, enhanced clearance headpiece

## SPECIFICATIONS

- Model Number: 1089
- Overall Length: 90"
- Overall Width (side rails up): 31.5"
- Weight Capacity: 500 lbs.
- Height Range (to litter top): 22.25" to 34"
- Patient Surface: 26" x 87"
- Caster Diameter: 8
- Trend/Reverse Trend.: ±18°

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)