

# Stryker 5051 Eye Surgery Stretcher Chair - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 5051 Eye Surgery Stretcher Chair is designed for surgery centers and space-limited environments, offering optimal patient positioning without compromising care quality. The chair features a dual-articulation, enhanced-clearance headpiece for surgical access, with easy adjustment from horizontal to upright positioning to support superior clinical performance and patient satisfaction.

This refurbished unit includes a new 2" pad set, 6" casters for smooth mobility, and dual-cylinder pneumatic backrest with dual pedestal hydraulics for reliable operation. The 24" wide patient surface accommodates various patient needs while maintaining a compact, user-friendly design.

The chair comes equipped with a thermoformed ABS hood featuring an integrated storage tray and oxygen bottle holder, along with swing-down siderails with comfort armrests and four IV receptacles for comprehensive surgical support.

## FEATURES

- Dual-Articulation, Enhanced-Clearance Headpiece
- 24" Wide Patient Surface
- Adjustable Height Footrest
- Dual Cylinder Pneumatic Backrest
- Dual Pedestal Hydraulics
- Thermoformed ABS Hood with Integrated Storage Tray
- Oxygen Bottle Holder
- Four-wheel, Steel-ring brakes with dual side activators

## SPECIFICATIONS

- Overall Length: 76"
- Overall Width: 30"
- Weight Capacity: 400 lbs.
- Height Range (to littler top): 22" to 33.5"
- Patient Surface: 24" x 74"
- Caster Diameter: 6"
- Litter Positioning: Backrest - 0° to 90°
- Trend./Reverse Trend. - ±18°

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)