

Stryker 1069 Eye Surgery Stretcher - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 1069 Eye Surgery Stretcher is engineered specifically for ophthalmic and head procedures, providing the greatest head-end clearance in the industry to give surgeons ample legroom for delicate operations. A firmly locking work surface and mobile base enhance both patient safety and clinical efficiency. The stretcher supports both top-down and side approaches with built-in controls in the head section for quick adjustments.

The dual-moving headpiece provides better access during procedures, while a flexible air tube manages airflow and holds drapes in place to maintain a sterile field. Dual-cylinder backrest support and dual pedestal hydraulics enable easy height adjustment. Side-mounted controls allow staff to adjust the stretcher from either side, making positioning simple and efficient.

Additional features include fold-down siderails, four IV receptacles, and four-wheel steel-ring brakes with foot-end brake and steer pedal for enhanced mobility and safety.

FEATURES

- New Patient Pad and Concave Headrest Pad
- Dual Cylinder Pneumatic Backrest
- Dual Pedestal Hydraulics
- Dual Side Hydraulic Controls
- Fold-down Siderails
- Four IV Receptacles
- Four-wheel Steel-ring Brakes with foot-end Brake & Steer Pedal
- Dual-articulation, Enhanced-Clearance Headpiece

SPECIFICATIONS

- Model Number: 1069
- Overall Length: 90"
- Overall Width (side rails up): 31.5"
- Weight Capacity: 500 lbs.
- Height Range (to litter top): 22.5" to 34"
- Patient Surface: 26" x 87"
- Headpiece: Height Adjustment - ± 4.5 ", Headpiece Tilt: ± 4.5 "
- Caster Diameter: 8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)