

# Stryker 1079 Eye Surgery Stretcher - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 1079 Eye Surgery Stretcher is designed specifically for ophthalmic procedures, offering exceptional head-end clearance that provides surgeons with extra legroom to work comfortably during delicate eye operations. The stretcher accommodates both superior and temporal surgical approaches with an adjustable headpiece featuring easy-to-use controls for precise patient positioning.

The unit features a secure locking surface to maintain patient stability and smooth-rolling design for effortless movement and repositioning in clinical settings. A flexible tube aids airflow while holding surgical drapes in place, and the adjustable backrest combined with dual-side hydraulic controls enables convenient operation from either side of the stretcher.

Built with a strong base that maintains steady height and fold-down siderails for accessibility, this refurbished stretcher is engineered for reliable performance in eye surgery applications. The unit includes four IV receptacles, four-wheel steel-ring brakes with foot-end brake and steer pedal, and integrated storage capabilities.

## FEATURES

- New Patient Pad and Concave Headrest Pad
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual-side hydraulic controls
- Fold-down side rails
- Four IV receptacles
- Four-wheel, steel-ring brakes with foot-end brake & steer pedal
- Dual-articulation, enhanced-clearance headpiece

## SPECIFICATIONS

- Model Number: 1079
- Overall Length: 90"
- Overall Width (side rails up): 31.5"
- Weight Capacity: 500 lbs.
- Height Range (to litter top): 23.25" to 36"
- Patient Surface: 26" x 87"
- Caster Diameter: 8"
- Headpiece: Height Adjustment -  $\pm 4.5"$ , Headpiece Tilt:  $\pm 45^\circ$

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)