

Stryker 1068 Eye Surgery Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker 1068 Eye Surgery Stretcher is specifically designed for ophthalmic procedures, offering precise patient positioning and optimal surgeon access. It features adjustable articulation including Trendelenburg/Reverse Trendelenburg at $\pm 18^\circ$, Fowler positioning from 0° to 90° , and a knee gatch for enhanced patient comfort and surgical positioning.

The stretcher provides a height range of 22" to 35.5", allowing easy and comfortable access for the surgeon. It is equipped with 8" casters for smooth maneuverability and 13" x 55" side rails for patient safety. The headrest offers $\pm 3^\circ$ height adjustment and $\pm 45^\circ$ tilt capability.

With a maximum patient weight capacity of 500 lbs and a patient surface width of 26", the Stryker 1068 combines durable construction with ergonomic features. Each unit comes with a new pad set including a concave headrest pad and includes a warranty.

FEATURES

- New Pad Set with Concave Headrest Pad
- Height Range: 22" to 35.5"
- 8" Casters
- Fowler: 0° to 90°
- Knee Gatch: 0° to 30°
- Head Piece Height Adjustment is $\pm 3^\circ$
- Trendelenburg/Head Elevations $\pm 18^\circ$
- Side Rails: 13" x 55"

SPECIFICATIONS

- Height Range: 22" to 35.5"
- 8" Diameter Casters
- Fowler: 0° to 90°
- Head Piece Height Adjustment is $\pm 3^\circ$
- Head Piece Tilt: $\pm 45^\circ$
- Trendelenburg/Head Elevations $\pm 18^\circ$
- Side Rails: 13" x 55"
- Maximum Weight Capacity: 500 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)