

Stryker 5050 Stretcher / Chair - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 5050 Stretcher/Chair, Refurbished. The refurbished Stryker 5050 Stretcher/Chair offers a compact, versatile design, ideal for medical facilities where space is limited. Its dual-purpose function allows it to convert easily from stretcher to chair, making it perfect for a range of clinical uses.

In addition, the CPR-capable frame and 400 lb. weight capacity provide reliability and support in various patient care scenarios.

Furthermore, the 5050 features a dual-cylinder pneumatic backrest, adjustable-height footrest, and dual pedestal hydraulics. Furthermore, additional features include hands-free Trendelenburg and Reverse Trendelenburg, a drop-down leg section, central brake and steer system, fixed-position push bar, four IV receptacles, fold-down side rails with comfort rest, and smooth-rolling 6" casters.

FEATURES

- Hands-free Trendelenburg/Reverse Trendelenburg
- Dual Cylinder Pneumatic Assist Back-Up/Down
- Height Up/Down
- Dual Hydraulic Jacks which optimizes stability for a more secure treatment platform
- Fold-Down Siderails

SPECIFICATIONS

- Overall Dimensions: 30"W x 76"L
- Height Range: 22" - 33.5"
- Back-Up/Down (Pneumatic Assist)
- Casters: 6"
- 400 lbs Weight Capacity

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)