

# SchureMed 11"L x 8"W Table Width Extender | 800-0081



Stock Photo - Options May Vary

## FACTORY NEW

The SchureMed Table Width Extender is a surgical positioning device designed to enhance patient safety and surgeon access during operating room procedures. It mounts securely on the O.R. table side rail and provides firm support during upper extremity and thoracic surgical procedures. The extension accommodates a 500 lb patient load and includes a deluxe pad for patient comfort.

The unique design allows any standard table accessory to be mounted to the extension, providing flexibility for various surgical setups. Height adjustment enables the extension to be leveled with the operating table, with or without an X-ray top. Installation is straightforward and can be completed by a single person.

The device is maintenance-free, latex-free, and built to last five years under normal use conditions. It weighs only 4 lbs, making it easy to handle while maintaining structural integrity for demanding surgical environments.

## FEATURES

- Surgical Table Bariatric Width Extension
- Mounts on the O.R. table side rail
- Provides firm support for the surgeon during upper extremity and thoracic surgical procedures
- Unique design allows any table accessory to be mounted to this extension
- Deluxe pad included
- Attaches to any O.R. side rail with mounting clamp (sold separately)
- Raise and lower the height of the extension to be level with the operating table, either with or without an X-ray top
- Single-person installation

## SPECIFICATIONS

- Length: 11" +/- 0.5"
- Width: 8" +/- 0.5"
- Depth: 8" +/- 0.5"
- Device Weight: 4 lbs +/- 0.5 lbs
- Weight Capacity: Device supports 500 lb proportional patient load (6'4" tall patient per 99% human body model)
- Product not made with (NRL) Natural Rubber Latex
- Product is maintenance-free
- The life of the device is five years under normal use

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)