

SchureMed 15"L X 8"W Table Width Extender | 800-0079



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed 15"L x 8"W Table Width Extender is a surgical table accessory designed to enhance patient positioning during operating room procedures. It mounts securely on the O.R. table side rail and provides firm support for the surgeon during upper extremity and thoracic surgical procedures. The extension features a unique design that allows any table accessory to be mounted to it, and includes a deluxe pad.

This bariatric table width extension raises and lowers to be level with the operating table, accommodating both standard and X-ray tops. It attaches to any O.R. side rail with a mounting clamp and requires a simple clamp for installation. The device is maintenance-free and supports a 500 lb proportional patient load.

The extension is designed for single-person installation and is not made with natural rubber latex. It carries a two-year warranty covering manufacturing defects, with an expected device life of five years under normal use.

FEATURES

- Surgical Table Bariatric Width Extension
- Mounts on the O.R. table side rail
- Provides firm support for the surgeon during upper extremity and thoracic surgical procedures
- Unique design allows any table accessory to be mounted to this extension
- Deluxe pad included
- Raise and lower the height of the extension to be level with the operating table; either with or without an X-ray top
- Product is maintenance-free
- The life of the device is five years under normal use

SPECIFICATIONS

- Device Dimensions (maximum): 15"L X 8" W
- Length: 15" +/- 0.5"
- Width: 8" +/- 0.5"
- Depth: 8" +/- 0.5"
- Device Weight: 4lbs +/- 0.5"
- Single-Person Installation
- The device supports a 500 lb proportional patient load (6'4" tall patient per 99% human body model)
- Product not made with (NRL) Natural Rubber Latex

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)