

SchureMed 20"L x 4" W Table Width Extender | 800-0077



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Table Width Extender is a surgical positioning device designed to enhance patient safety and surgeon access during operating room procedures. It mounts securely to the O.R. table side rail and provides firm support during upper extremity and thoracic surgical procedures. The extender accommodates any standard table accessory and includes a deluxe pad.

The extension height adjusts to align with the operating table surface, whether or not an X-ray top is in use. Its unique design allows flexible mounting of additional accessories to meet specific procedural requirements. Installation is straightforward and can be completed by a single person.

The device attaches to any O.R. side rail using a mounting clamp (sold separately). This positioning solution helps optimize surgical access and patient comfort during a variety of procedures.

FEATURES

- Surgical Table Bariatric Width Extension Mounts on the O.R. table side rail
- Provides firm support for the surgeon during upper extremity and thoracic surgical procedures
- Unique design allows any table accessory to be mounted to this extension
- Deluxe pad included
- Attaches to any O.R. side rail with mounting clamp (sold separately)
- Raise and lower the height of the extension to be level with the operating table; either with or without an X-ray top
- Required Simple Clamp (sold separately, see Option Tab)

SPECIFICATIONS

- Dimensions: 20" L X 4" W
- Length: 20" +/- 0.5"
- Width: 4" +/- 0.5"
- Depth: 8" +/- .05"
- Device Weight: 4.5 +/- 0.5"
- Single-person installation

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)