

# SchureMed Schure Socket XL | 800-0134



Stock Photo - Options May Vary

## FACTORY NEW

The Schure Socket XL is a universal mounting socket designed for secure attachment of surgical accessories to operating room tables. It snaps onto any standard OR table side rail and accommodates round mounting posts up to 5/8" diameter as well as rectangular posts up to 1" x 3/8". The socket features a hardened clamping surface that creates maximum clamping force to prevent movement of connected accessories.

The quick-turn design allows secure connection and disconnection at any point along the side rail, including over side rail standoffs. Installation is straightforward and can be completed by a single person. The socket is compatible with most OR table designs and brands, making it a versatile solution for affixing a wide range of surgical table accessories.

The device supports a patient weight capacity of 800 lbs and is built to withstand the demands of gynecological, urological, laparoscopic, and orthopedic surgical procedures. The socket is maintenance-free and carries a two-year warranty against manufacturing defects.

## FEATURES

- Clark Sockets by SchureMed are designed for a wide range of purposes; affixing most operating room table accessories to most OR table designs and brands
- Quickly and easily connects and disconnects securely with a quick turn anywhere along the side rail, even over side rail standoff
- The hardened clamping surface creates maximum clamping force stopping all movement of operating room accessories
- Accepts any round mounting posts up to 5/8" D (1.6 cm) and rectangular mounting posts up to 1" x 3/8" (2.5 cm x .95 cm)

## SPECIFICATIONS

- Compatible with all standard operating room table side rails (3/8"W x 1-1/8"H)
- Length: 4" +/- 0.5"
- Width: 2.25" +/- 0.5"
- Depth: 2.25" +/- 0.5"
- Device Weight: 1.5 +/- 0.5 lbs
- Patient weight capacity is 800 lbs
- The product warranty covers the clamp from manufacturing defects for two years
- The life of the device is five years under normal use

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)