

SchureMed Elite 40° Trendelenburg System | 804-0118



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Elite 40° Trendelenburg System (804-0118) addresses key challenges in the Trendelenburg position, including patient slippage, nerve injury, and decubitus sore formation. The system features adjustable angles for precise positioning to accommodate various surgical needs and innovative padded, non-slip surfaces that securely hold the patient, reducing shifting and injury risk during surgery.

The Elite Trendelenburg System is compatible with most operating tables and integrates easily into existing surgical setups. Its durable construction and ease of use make it a reliable tool for optimizing patient care in high-risk procedures, ultimately enhancing surgical outcomes and reducing complications.

FEATURES

- Exceptionally comfortable gel and foam combination, which offers restraint-free pressure relief and stabilization of the patient.
- A unique dual-layer, integrated gel structure conforms to the patient's shape, improving immersion and comfort.
- A combination of memory foam and latticed gel cells provides maximum pressure relief.
- Advanced carbon composite shoulder supports provide additional comfort and stabilization for steeper procedures.
- Surgical Positions: Recommended for robotic-assisted surgery, abdominal surgery & gynecological surgery.

SPECIFICATIONS

- 12-Month Warranty



**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com