

# SchureMed Foot Extension | 800-0036



Stock Photo - Options May Vary

## FACTORY NEW

The SchureMed Foot Extension (800-0036) is a versatile surgical table accessory designed to enhance patient positioning during operating room procedures. It functions both horizontally as a table extension and vertically as foot support for reverse Trendelenburg positions, providing flexibility for optimal surgical access and patient comfort.

The extension is completely radiolucent, allowing unobstructed imaging during procedures, and fits all standard O.R. tables. A deluxe 2-inch thick pad is included for patient comfort and protection. The device is latex-free and requires no maintenance during its five-year service life.

With a 500-pound patient load capacity, the foot extension supports a wide range of patient sizes. Installation is straightforward and can be completed by a single person. Two simple clamps (sold separately) are required for attachment to the surgical table rail.

## FEATURES

- Can be used horizontally as a table extension
- Can be used vertically as foot support for reverse Trendelenburg positions Deluxe pad included.
- Fits all O.R. tables
- Completely Radiolucent
- Product not made with (NRL) Natural Rubber Latex
- The device supports a 500 lb proportional patient load (6'4" tall patient per 99% human body model)
- Product is maintenance-free
- The life of the device is five years under normal use

## SPECIFICATIONS

- Overall Dimensions: 20"W X 10"L
- Pad: 2" thick, included
- Completely radiolucent
- Length: 18.5" +/- 0.5"
- Width: 24" +/- 0.5" (Extended), 21.750" +/- 0.5" (Collapsed)
- Depth: 3.75" +/- 0.5" (with pad)
- Device Weight: 5.25 +/- 0.5 lbs
- Single-person installation

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)