

SchureMed Hand Shield | 807-0010



Stock Photo - Options May Vary

FACTORY NEW

SchureMed Hand Shield: 807-0010. Enjoy full topside protection without sacrificing ease of movement.

The lighter-weight Hand Shield's open palm design allows free movement of fingers and thumbs, making it ideal for procedures that require dexterity. It's especially useful for holding the patient during radiographic procedures, ensuring both safety and comfort.

The ergonomic design provides a secure fit while minimizing fatigue, enabling longer, more efficient use. This durable, high-quality shield offers essential protection while maintaining flexibility and ease of use in critical surgical settings.

FEATURES

- Thumb strap for improved control
- Velcro strap closure for quick on and off
- Manufactured with soft lead for superior flexibility
- Sold as a pair (set of 2)

SPECIFICATIONS

- Closure: Velcro
- Provides 0.5mm Pb Equiv protection

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)