

# SchureMed Head Lead | Reusable



*Stock Photo - Options May Vary*

## **FACTORY NEW**

SchureMed Head Lead: Reusable. Put a lid on radiation exposure with a lightweight, protective, and reusable cap. Reusable Head Lead is made with non-leaded Earth-Lite™ or leaded Ultra-Lite™ material and is available in two different levels of protection: .125 mm Pb & .25 mm Pb lead equivalency.

The inside is constructed in a wick-away fabric that will keep you cool & comfortable while wearing.

Three sizes are available; all can be adjusted with the rear drawstring.

## **FEATURES**

- Lightweight, non-leaded material
- Moisture-wicking fabric
- Adjustable elastic drawstring
- Three sizes to fit all heads
- Three sizes available: Small (21" circumference x 4" deep),

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)