

# SchureMed Lateral Patient Positioners | 800-0043



Stock Photo - Options May Vary

## FACTORY NEW

SchureMed Lateral Patient Positioners (800-0043) are designed to address critical patient positioning challenges in gynecological, urological, laparoscopic, and orthopedic surgical procedures. These positioners clamp quickly to standard operating room tables and feature three swivel braces that conform to any size patient while maintaining complete stability.

The system is simple to set up and offers an effective alternative for lateral positioning during surgery. Three swivel positioners with deluxe foam pads rotate to accommodate various patient sizes, providing ease of use and reliable support throughout the procedure.

The positioners feature a 5/8" stainless steel mounting post and are sold in sets. Three Schure Socket XLs are required for operation and are sold separately.

## FEATURES

- Simple and effective lateral patient positioners
- Clamp quickly to standard operating room table
- Three swivel positioners conform to any size patient
- Sold in sets
- Simple to set up and easy-to-use
- 5/8" Stainless Steel Mounting Post
- Deluxe foam pads included
- Three required Schure Socket XLs are sold separately

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)