

SchureMed Lead Gloves | 807-0064



Stock Photo - Options May Vary

FACTORY NEW

SchureMed Lead Gloves: 807-0064. You don't have to sacrifice comfort or dexterity with these leaded finger gloves. Seamless construction provides a functional design that delivers a better fit and flexibility.

The gloves offer superior protection against radiation while allowing for natural hand movement. Their lightweight material reduces hand fatigue, making them ideal for extended use in radiographic procedures. The ergonomic fit ensures maximum comfort, while the durable design withstands daily wear and tear.

Perfect for professionals who need both safety and dexterity during high-precision tasks.

FEATURES

- Seamless design
- Vinyl-dipped
- Brushed nylon interior
- One size fits all

SPECIFICATIONS

- One size fits all
- 0.5 mm lead equivalent (standard configuration provides 97% direct beam attenuation at 100 kVp)
- 15" overall length

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)