

SchureMed Narrow Head Rest | 800-0052



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Narrow Head Rest (800-0052) is an operating room table extension that replaces the standard headrest to improve surgical access during head and neck procedures. The unique tapered design narrows from 20 inches wide to 8 inches at the head, providing optimal visibility and workspace for surgeons performing ENT, plastic surgery, neurosurgery, or ophthalmic procedures. A deluxe pressure management pad is included with the headrest.

FEATURES

- Fits all O.R. tables
- 20"W reduces to 8"W at the head.
- The radiolucent headrest includes a 2" premium pressure management pad
- Wrist Rest adjusts 14" vertically and the inside diameter width adjusts 5" to 11" (not included, see Options)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)