

SchureMed Oasis Elite Prone Head Rest

| 804-0025



Stock Photo - Options May Vary

FACTORY NEW

The Oasis Elite Prone Head Rest is anatomically designed for comfortable face-down positioning with respiratory tubes in place during prone surgery. It provides excellent support, stability, comfort, and protection with openings on both sides to accommodate breathing apparatus. The product is lightweight, durable, and easy to handle for operating room staff.

FEATURES

- Anatomically designed so that the patient's head can comfortably rest face down with respiratory tubes in place
- Provides support, stability, comfort, and protection during surgery
- In the prone position, with openings on both sides
- Lightweight, durable, and easy to handle

SPECIFICATIONS

- Adult Size Dimensions: 9.45"W x 11.02"L x 5.51"H
- Adolescent/Adult Dimensions: 8.66"W x 10.75"L x 5.31"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)