

SchureMed Oasis Elite Prone Positioner

| 804-0039

FACTORY NEW

The SchureMed Oasis Elite Prone Positioner is designed to support patients during surgical procedures performed in the prone position. It features a gel and foam combination that provides effective positioning while minimizing pressure points to enhance patient safety and comfort. The durable design allows for easy integration into various surgical settings, making it an essential tool for prone positioning procedures.



Stock Photo - Options May Vary

SPECIFICATIONS

- Dimensions: 20.87"W x 24.02"L x 5.12"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)