

# SchureMed Patient Positioner Lateral Braces | 800-0042



Stock Photo - Options May Vary

## FACTORY NEW

SchureMed Lateral Braces are surgical patient positioning devices designed to provide total lateral stability during gynecological, urological, laparoscopic, and orthopedic procedures. The braces clamp quickly to standard operating room tables and conform to any size patient, making them suitable for total hip replacements and other procedures requiring secure lateral positioning.

These positioning braces feature a 5/8" stainless steel mounting post and are simple to set up and use in the operating room. Three multi-axis braces work together to easily accommodate patients of varying sizes while maintaining perfect lateral alignment.

The system includes deluxe foam pads with protective covers that provide extreme pressure relief and patient comfort during extended surgical procedures. The braces are sold in sets and require three Schure Socket XLs for mounting, which are sold separately.

## FEATURES

- Sold in sets
- Simple to set up and easy-to-use
- Clamps quickly to all standard U.S. operating tables
- 5/8" Stainless Steel Mounting Post
- Secure and comfortable positioning for total hip replacement, and any other procedure requiring total lateral stability
- Three multi-axis braces easily accommodate any size patient for perfect lateral positioning
- Pressure relief foam with Deluxe cover provide extreme protection for the patient
- Deluxe foam pads included

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)