

# SchureMed Powered Beach Chair | 800-0004



Stock Photo - Options May Vary

## FACTORY NEW

The SchureMed Powered Beach Chair (800-0004) is a shoulder surgery positioning system featuring a self-contained motor independent of the surgical table. The double-ball joint headrest allows the patient's head to tilt 30° in all directions—forward, backward, and side-to-side—accommodating patients with kyphosis and ensuring precise positioning for optimal surgical outcomes.

The chair adjusts from 0-90° in 8 seconds with controlled, motorized movement via push-button motion-control technology. It supports arthroscopy patients up to 500 lbs with a total lift capacity of 1,350 lbs. The unit attaches to any OR table using adjustable mounting clamps and includes a lateral brace, shoulder pull-away panels for surgical exposure, and a five-piece deluxe pad set for patient comfort.

The lateral brace attaches directly to the beach chair rather than the OR table, ensuring safer patient positioning. The system eliminates manual positioning requirements and provides full 0-90° articulation range for flexible surgical access.

## FEATURES

- The self-contained motor that does not rely on the surgical table for power
- A double-ball joint headrest allows the patient's head to tilt 30° in every direction, forward, backward, and side-to-side
- Adjusts from 0-90° in 8 seconds in a slow, controlled movement
- Motion-control technology eliminates manual patient positioning with a simple push-button hand control
- 1,350 lbs. of lift securely supports arthroscopy patients up to 500 lbs.
- Attaches to ANY operating room table with our integrated, adjustable mounting clamps
- Lateral brace attaches to the beach chair and not O.R. table side rail
- Shoulder pull-away panels provide excellent surgical exposure

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)